## CORMAY AMMONIA / ETHANOL CONTROLS

### CORMAY AMMONIA / ETHANOL CONTROLS

Nr kat. 5-163 3 x 3 ml

- AMMONIA / ETHANOL CONTROL 1 1 x 3 ml
- AMMONIA / ETHANOL CONTROL 2 1 x 3 ml
- AMMONIA / ETHANOL CONTROL 3 1 x 3 ml

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#### INTENDED USE

CORMAY AMMONIA / ETHANOL CONTROLS kit is intended for quality control of ethanol and ammonia determinations.

Controls are prepared from bovine serum albumin with preservatives and stabilizers added.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

#### STABILITY AND STORAGE

- Unopened controls are stable at 2-8°C until the expiration date printed on the package.
- Stability after opening: 20 days at 2-8°C tightly capped.

#### PROCEDURE

Controls are ready to use.

- Allow the controls to equilibrate to room temperature prior to use.
- Mix controls gently and accurately before use.
- Do not allow controls to stand for long time periods at room temperature. After use store vials tightly closed in refrigerator.

#### WARNINGS AND NOTES

- CORMAY AMMONIA / ETHANOL CONTROLS meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

CORMAY AMMONIA / ETHANOL CONTROLS contain reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one, mixture (3:1).

Warning:

H317 May cause an allergic skin reaction.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice.

P363 Wash contaminated clothing before reuse.

- Do not use controls after expiry date
- Do not interchange the caps and plugs.
- Turbidity of controls may signify a contamination. Do not use turbid controls.
- Controls are not intended for use as a standard.

### ASSIGNED VALUES

AMMONIA / ETHANOL CONTROL 1			Lot: 807-41
COMPONENT	MEAN		RANGE
Ammonia	61,3	µg/dl	42,6 – 80,1
	36,0	µmol/l	25,0 – 47,0
Ethyl alcohol	42,0	mg/dl	34,0 – 50,0
	9,11	mmol/l	7,38 – 10,9
AMMONIA / ETHANOL CONTROL 2			Lot: 807-41
COMPONENT	MEAN		RANGE
Ammonia	186	µg/dl	148 – 223
	109	µmol/l	87,0 – 131
Ethyl alcohol	103	mg/dl	82,0 – 124
	22,4	mmol/l	17,8 – 26,9
AMMONIA / ETHANOL CONTROL 3			Lot: 807-43
COMPONENT	MEAN		RANGE
Ammonia	489	µg/dl	392 – 586
	287	µmol/l	230 – 344
Ethyl alcohol	250	mg/dl	200 – 300
	54,3	mmol/l	43,4 – 65,1

Assigned values are specific for lot number of control material  
 It is recommended that each laboratory establish its own means and acceptable ranges.

#### WASTE MANAGEMENT

Please refer to local legal requirements

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