CORMAY D-DIMER CALIBRATOR

(PL)

CORMAY D-DIMER CALIBRATOR

Nr kat. 4-259 4 x 1 ml

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Cat. No 4-259 4 x 1 ml

INTENDED USE

CORMAY D-DIMER CALIBRATOR is designed for use as a calibrator for determination of D-dimer concentration. The exact concentration of D-dimer is printed on the label and in the table below.

The calibrators must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

FDPs (Fibrin and Fibrinogen Degradation Products) are generic name for several degradation products that are formed as the result of plasmin mediated, proteolytic degradation of fibrin and fibrinogen. FDP assay become important test to diagnose or monitor fibrinolytic disorder, especially disseminated intravascular coagulation (DIC). But FDP assay also detects the degradation products of fibrinogenolysis. D-dimer assay is more specific for fibrinolysis because it can only be produced as the result of plasmin mediated, proteolytic degradation of fibrin.

RECONSTITUTION AND PREPARATION

Carefully remove the rubber cap from vial and add 1 ml D-DIMER DILUENT to the vial. Recap, swirl gently and let stand undisturbed for 15 minutes at room temperature.

Afterwards, using D-DIMER DILUENT, prepare a series of dilutions, according to the table below:

CORMAY D-DIMER CALIBRATOR						Lot: 807-45
Dilution	1/32	1/16	1/8	1/4	1/2	neat
D-dimer concentration (µg/ml FEU)	0.65	1.3	2.6	5.2	10.4	20.8

In case of analysis on **Hitachi 911/912** analyser, dilution of calibrator are made automatically. Calibrators values for determinations on Hitachi 911/912:

CORMAY D-DIMER CALIBRATOR		Lot: 807-45
Calibrator	Concentration (µg/ml FEU)	
1	0.65	
2	1.3	
3	2.6	
4	5.2	
5	10.4	
6	20.8	

STABILITY AND STORAGE

1. Unopened calibrator is stable at 2-10°C until the expiration date printed on the vial.
2. Once reconstituted, calibrator is stable for 2 weeks when stored at 2-10°C.

WARNINGS AND NOTES

1. This calibrator has been tested for the HIV antibody, HBsAg and HCV and found to be non-reactive.

However this material should be handled as though capable of transmitting infectious disease.

2. D-DIMER DILUENT is present in kits for determination of D-dimer concentration.
3. Do not mix different lot.
4. Do not use expired lot.
5. D-DIMER DILUENT contains < 0.1% sodium azide as a preservative. Avoid contact with skin and mucous membranes.

WASTE MANAGEMENT

Please refer to local legal requirements.

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