## CORMAY HbA<sub>1c</sub> DIRECT CALIBRATORS

### CORMAY HbA<sub>1c</sub> DIRECT CALIBRATORS (PL)

- Nr kat. 4-308 4 x 0,5 ml
- HbA<sub>1c</sub> DIRECT CALIBRATOR 1 1 x 0,5 ml
  - HbA<sub>1c</sub> DIRECT CALIBRATOR 2 1 x 0,5 ml
  - HbA<sub>1c</sub> DIRECT CALIBRATOR 3 1 x 0,5 ml
  - HbA<sub>1c</sub> DIRECT CALIBRATOR 4 1 x 0,5 ml

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#### INTENDED USE

CORMAY HbA<sub>1c</sub> DIRECT CALIBRATORS kit is dedicated for preparation of the calibration curve for direct method HbA<sub>1c</sub> determinations.

The calibrators must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

#### INTRODUCTION

The determination of HbA<sub>1c</sub> is most commonly performed for the evaluation of glycemic control in diabetes mellitus. HbA<sub>1c</sub> values provide an indication of glucose levels over the preceding 4-8 weeks.

#### RECONSTITUTION

Carefully remove the cap from vial and add exactly 0.5 ml of deionized water. Recap, gently mix for 10 minutes until all material has dissolved.

#### STABILITY AND STORAGE

1. Unopened calibrator is stable at 2-8°C until the expiration date printed on the vial. Protect from light and heat!
2. The constituent in reconstituted calibrator is stable for: 4 weeks at 2-8°C or 8 weeks at -20°C.

#### WARNINGS AND NOTES

- The calibrator is based on human blood that has been tested for HBsAg and antibodies to HIV and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
- Avoid contact with skin and mucous membranes.
- Before the assay HbA<sub>1c</sub> DIRECT reconstituted calibrators should be treated with HEMOLYSING REAGENT according to HbA<sub>1c</sub> DIRECT kit insert.

### ASSIGNED VALUES

Values in NGSP units [%]:

CORMAY HbA <sub>1c</sub> DIRECT CALIBRATORS		Lot: 806-26			
Analyser	Calibrator				
	1	2	3	4	
Prestige 24i Bolis 24i Premium Bolis 15i Bolis 50i ACCENT-200 ACCENT-200 II GEN ACCENT-220S	5.9%	8.5%	11.3%	15.4%	
BS-400, BS-480 Olympus AU400/AU640 BS-800, BS-800M	5.2%	8.1%	11.1%	15.6%	
ACCENT-300 HITACHI 911/912	5.3%	8.6%	11.4%	16.0%	

#### TRACEABILITY

HbA<sub>1c</sub> DIRECT calibrators are traceable to IFCC standards through the Tosoh calibrator.

#### WASTE MANAGEMENT

Please refer to local legal requirements.

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