## CORMAY AMMONIA / ETHANOL CONTROLS

|                                   |          |
|-----------------------------------|----------|
| Nr kat. 5-163                     | (PL)     |
| CORMAY AMMONIA / ETHANOL CONTROLS | 3 x 3 ml |
| AMMONIA / ETHANOL CONTROL 1       | 1 x 3 ml |
| AMMONIA / ETHANOL CONTROL 2       | 1 x 3 ml |
| AMMONIA / ETHANOL CONTROL 3       | 1 x 3 ml |

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### INTENDED USE

CORMAY AMMONIA / ETHANOL CONTROLS kit is intended for quality control of ethanol and ammonia determinations.

Controls are prepared from bovine serum albumin with chemicals, preservatives and stabilizers added.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### STABILITY AND STORAGE

- Unopened controls are stable at 2-8°C until the expiration date printed on the package.
- Stability after opening: 20 days at 2-8°C, tightly capped.

### PROCEDURE

Controls are ready to use.

- Allow the controls to equilibrate to room temperature prior to use!
- Mix controls gently and accurately before use.
- Do not allow controls to stand for long time periods at room temperature. After use store vials tightly closed in refrigerator.

### WARNINGS AND NOTES

- Do not use controls after expiry date
- Do not interchange the caps and plugs.
- Turbidity of controls may signify a contamination. Do not use turbid controls.
- Controls are not intended for use as a standard.
- CORMAY AMMONIA / ETHANOL CONTROLS meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

CORMAY AMMONIA / ETHANOL CONTROLS contains 5-chloro-2-methyl-2H-isothiazol-3-one.

#### Warning:

H317 May cause an allergic skin reaction.

P280 Wear protective gloves.

P261: Avoid breathing spray.

P363 Wash contaminated clothing before reuse.

P302+P352 IF ON SKIN: Wash with plenty

of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice.

### ASSIGNED VALUES

| AMMONIA / ETHANOL CONTROL 1 |      |        | Lot: 305-60 |
|-----------------------------|------|--------|-------------|
| COMPONENT                   | MEAN |        | RANGE       |
| Ammonia                     | 85.2 | µg/dl  | 59.6 – 111  |
|                             | 50.0 | µmol/l | 35.0 – 65.0 |
| Ethyl alcohol               | 28.0 | mg/dl  | 22.0 – 34.0 |
|                             | 6.08 | mmol/l | 4.77 – 7.38 |
| AMMONIA / ETHANOL CONTROL 2 |      |        | Lot: 305-60 |
| COMPONENT                   | MEAN |        | RANGE       |
| Ammonia                     | 210  | µg/dl  | 167 – 252   |
|                             | 123  | µmol/l | 98.0 – 148  |
| Ethyl alcohol               | 83.0 | mg/dl  | 66.0 – 100  |
|                             | 18.0 | mmol/l | 14.3 – 21.7 |
| AMMONIA / ETHANOL CONTROL 3 |      |        | Lot: 305-63 |
| COMPONENT                   | MEAN |        | RANGE       |
| Ammonia                     | 501  | µg/dl  | 400 – 601   |
|                             | 294  | µmol/l | 235 – 353   |
| Ethyl alcohol               | 237  | mg/dl  | 190 – 284   |
|                             | 51.4 | mmol/l | 41.2 – 61.6 |

Assigned values are specific for lot number of control material. It is recommended that each laboratory establish its own means and acceptable ranges.

### WASTE MANAGEMENT

Please refer to local legal requirements

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