

CORMAY IMMUNO-CONTROL II (L)

Nr kat. / Cat. No 4-290

Seria / Lot: 509-42

Data ważn. / Exp.: 2016.12

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
<b>BETA 2-MIKROGLOBULINA</b> <b>BETA 2-MICROGLOBULIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	1.70 mg/l	1.45 – 1.96 mg/l
<b>MIOGLOBINA</b> <b>MYOGLOBIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	103 ng/ml	88 – 118 ng/ml
<b>FERRYTYNA</b> <b>FERRITIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	101 ng/ml	86 – 116 ng/ml
<b>CAŁKOWITE IgE</b> <b>TOTAL IgE</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	66 IU/ml	50 – 83 IU/ml

CORMAY IMMUNO-CONTROL II (L)

Nr kat. / Cat. No 4-290

Seria / Lot: 509-42

Data ważn. / Exp.: 2016.12

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
<b>BETA 2-MIKROGLOBULINA</b> <b>BETA 2-MICROGLOBULIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	1.70 mg/l	1.45 – 1.96 mg/l
<b>MIOGLOBINA</b> <b>MYOGLOBIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	103 ng/ml	88 – 118 ng/ml
<b>FERRYTYNA</b> <b>FERRITIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	101 ng/ml	86 – 116 ng/ml
<b>CAŁKOWITE IgE</b> <b>TOTAL IgE</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	66 IU/ml	50 – 83 IU/ml

CORMAY IMMUNO-CONTROL II (H)

Nr kat. / Cat. No 4-290

Seria / Lot: 509-43

Data ważn. / Exp.: 2016.12

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
<b>BETA 2-MIKROGLOBULINA</b> <b>BETA 2-MICROGLOBULIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	7.73 mg/l	6.96 – 8.50 mg/l
<b>MIOGLOBINA</b> <b>MYOGLOBIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	299 ng/ml	269 – 329 ng/ml
<b>FERRYTYNA</b> <b>FERRITIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	394 ng/ml	335 – 453 ng/ml
<b>CAŁKOWITE IgE</b> <b>TOTAL IgE</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	406 IU/ml	345 – 467 IU/ml

CORMAY IMMUNO-CONTROL II (H)

Nr kat. / Cat. No 4-290

Seria / Lot: 509-43

Data ważn. / Exp.: 2016.12

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
<b>BETA 2-MIKROGLOBULINA</b> <b>BETA 2-MICROGLOBULIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	7.73 mg/l	6.96 – 8.50 mg/l
<b>MIOGLOBINA</b> <b>MYOGLOBIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	299 ng/ml	269 – 329 ng/ml
<b>FERRYTYNA</b> <b>FERRITIN</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	394 ng/ml	335 – 453 ng/ml
<b>CAŁKOWITE IgE</b> <b>TOTAL IgE</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	406 IU/ml	345 – 467 IU/ml