

CORMAY IMMUNO-CONTROL II (L)

Nr kat. / Cat. No 4-290

Seria / Lot: 612-45

Data ważn. / Exp.: 2018.06

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BETA 2-MIKROGLOBULINA BETA 2-MICROGLOBULIN	TIA-wzmacniane lateksem Latex-enhanced TIA	1.68 mg/l	1.43 – 1.93 mg/l
MIOGLOBINA MYOGLOBIN	TIA-wzmacniane lateksem Latex-enhanced TIA	102 ng/ml	87 – 117 ng/ml
FERRYTYNA FERRITIN	TIA-wzmacniane lateksem Latex-enhanced TIA	100 ng/ml	85 – 115 ng/ml
CAŁKOWITE IgE TOTAL IgE	TIA-wzmacniane lateksem Latex-enhanced TIA	65 IU/ml	49 – 81 IU/ml

CORMAY IMMUNO-CONTROL II (L)

Nr kat. / Cat. No 4-290

Seria / Lot: 612-45

Data ważn. / Exp.: 2018.06

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BETA 2-MIKROGLOBULINA BETA 2-MICROGLOBULIN	TIA-wzmacniane lateksem Latex-enhanced TIA	1.68 mg/l	1.43 – 1.93 mg/l
MIOGLOBINA MYOGLOBIN	TIA-wzmacniane lateksem Latex-enhanced TIA	102 ng/ml	87 – 117 ng/ml
FERRYTYNA FERRITIN	TIA-wzmacniane lateksem Latex-enhanced TIA	100 ng/ml	85 – 115 ng/ml
CAŁKOWITE IgE TOTAL IgE	TIA-wzmacniane lateksem Latex-enhanced TIA	65 IU/ml	49 – 81 IU/ml

CORMAY IMMUNO-CONTROL II (H)

Nr kat. / Cat. No 4-290

Seria / Lot: 612-46

Data ważn. / Exp.: 2018.06

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BETA 2-MIKROGLOBULINA BETA 2-MICROGLOBULIN	TIA-wzmacniane lateksem Latex-enhanced TIA	7.81 mg/l	7.03 – 8.59 mg/l
MIOGLOBINA MYOGLOBIN	TIA-wzmacniane lateksem Latex-enhanced TIA	305 ng/ml	275 – 336 ng/ml
FERRYTYNA FERRITIN	TIA-wzmacniane lateksem Latex-enhanced TIA	385 ng/ml	327 – 443 ng/ml
CAŁKOWITE IgE TOTAL IgE	TIA-wzmacniane lateksem Latex-enhanced TIA	404 IU/ml	343 – 465 IU/ml

CORMAY IMMUNO-CONTROL II (H)

Nr kat. / Cat. No 4-290

Seria / Lot: 612-46

Data ważn. / Exp.: 2018.06

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BETA 2-MIKROGLOBULINA BETA 2-MICROGLOBULIN	TIA-wzmacniane lateksem Latex-enhanced TIA	7.81 mg/l	7.03 – 8.59 mg/l
MIOGLOBINA MYOGLOBIN	TIA-wzmacniane lateksem Latex-enhanced TIA	305 ng/ml	275 – 336 ng/ml
FERRYTYNA FERRITIN	TIA-wzmacniane lateksem Latex-enhanced TIA	385 ng/ml	327 – 443 ng/ml
CAŁKOWITE IgE TOTAL IgE	TIA-wzmacniane lateksem Latex-enhanced TIA	404 IU/ml	343 – 465 IU/ml