

CORMAY IMMUNO-CONTROL II (L)

Nr kat. / Cat. No 4-290

Seria / Lot: 804-50

Data ważn. / Exp.: 2019.06

| PARAMETR / COMPONENT | METODA / METHOD | WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE | ZAKRES / RANGE |
|---|---|-----------------------------------|------------------|
| BETA 2-MIKROGLOBULINA BETA 2-MICROGLOBULIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 1.81 mg/l | 1.54 – 2.08 mg/l |
| MIOGLOBINA MYOGLOBIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 97 ng/ml | 82 – 112 ng/ml |
| FERRYTYNA FERRITIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 103 ng/ml | 88 – 118 ng/ml |
| CAŁKOWITE IgE TOTAL IgE | TIA-wzmacniane lateksem Latex-enhanced TIA | 67 IU/ml | 50 – 84 IU/ml |

CORMAY IMMUNO-CONTROL II (L)

Nr kat. / Cat. No 4-290

Seria / Lot: 804-50

Data ważn. / Exp.: 2019.06

| PARAMETR / COMPONENT | METODA / METHOD | WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE | ZAKRES / RANGE |
|---|---|-----------------------------------|------------------|
| BETA 2-MIKROGLOBULINA BETA 2-MICROGLOBULIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 1.81 mg/l | 1.54 – 2.08 mg/l |
| MIOGLOBINA MYOGLOBIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 97 ng/ml | 82 – 112 ng/ml |
| FERRYTYNA FERRITIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 103 ng/ml | 88 – 118 ng/ml |
| CAŁKOWITE IgE TOTAL IgE | TIA-wzmacniane lateksem Latex-enhanced TIA | 67 IU/ml | 50 – 84 IU/ml |

CORMAY IMMUNO-CONTROL II (H)

Nr kat. / Cat. No 4-290

Seria / Lot: 804-51

Data ważn. / Exp.: 2019.06

| PARAMETR / COMPONENT | METODA / METHOD | WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE | ZAKRES / RANGE |
|---|---|-----------------------------------|------------------|
| BETA 2-MIKROGLOBULINA BETA 2-MICROGLOBULIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 7.84 mg/l | 7.06 – 8.62 mg/l |
| MIOGLOBINA MYOGLOBIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 299 ng/ml | 269 – 329 ng/ml |
| FERRYTYNA FERRITIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 403 ng/ml | 343 – 463 ng/ml |
| CAŁKOWITE IgE TOTAL IgE | TIA-wzmacniane lateksem Latex-enhanced TIA | 409 IU/ml | 348 – 470 IU/ml |

CORMAY IMMUNO-CONTROL II (H)

Nr kat. / Cat. No 4-290

Seria / Lot: 804-51

Data ważn. / Exp.: 2019.06

| PARAMETR / COMPONENT | METODA / METHOD | WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE | ZAKRES / RANGE |
|---|---|-----------------------------------|------------------|
| BETA 2-MIKROGLOBULINA BETA 2-MICROGLOBULIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 7.84 mg/l | 7.06 – 8.62 mg/l |
| MIOGLOBINA MYOGLOBIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 299 ng/ml | 269 – 329 ng/ml |
| FERRYTYNA FERRITIN | TIA-wzmacniane lateksem Latex-enhanced TIA | 403 ng/ml | 343 – 463 ng/ml |
| CAŁKOWITE IgE TOTAL IgE | TIA-wzmacniane lateksem Latex-enhanced TIA | 409 IU/ml | 348 – 470 IU/ml |