

CORMAY AMMONIA / ETHANOL CONTROLS



KARTA WARTOŚCI WYZNACZONYCH / ASSAY SHEET

F07-QAD-01_05 15.12.2023

Nr kat. / Cat. No	Seria / Lot	Data ważności / Expiration date
5-163	403-636366	2025-03

PARAMETR / PARAMETER	KONTROLA / CONTROL	WARTOŚĆ ŚRODKOWA / TARGET VALUE	ZAKRES / RANGE
AMONIAK / AMMONIA	AMMONIA / ETHANOL CONTROL 1	86.9 µg/dL	61.3 - 112 µg/dL
		51.0 µmol/L	36.0 – 66.0 µmol/L
	AMMONIA / ETHANOL CONTROL 2	204 µg/dL	164 - 245 µg/dL
		120 µmol/L	96.0 - 144 µmol/L
	AMMONIA / ETHANOL CONTROL 3	477 µg/dL	382 - 572 µg/dL
		280 µmol/L	224 - 336 µmol/L
ALKOHOL ETYLOWY / ETHYL ALCOHOL	AMMONIA / ETHANOL CONTROL 1	38.0 mg/dL	30.0 – 46.0 mg/dL
		8.25 mmol/L	6.51 – 9.98 mmol/L
	AMMONIA / ETHANOL CONTROL 2	94.0 mg/dL	75.0 - 113 mg/dL
		20.4 mmol/L	16.3 – 24.5 mmol/L
	AMMONIA / ETHANOL CONTROL 3	242 mg/dL	194 - 290 mg/dL
		52.5 mmol/L	42.1 – 62.9 mmol/L