

CORMAY IMMUNO-CONTROL I (L)

Nr kat. / Cat. No 4-288

Seria / Lot: 210-72

Data ważn. / Exp.: 2024.07

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
<b>BIALKO C-REAKTYWNE (CRP)</b> <b>C-REACTIVE PROTEIN (CRP)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	0.90 mg/dl	0.63 – 1.17 mg/dl
<b>ANTY-STREPTOLIZYNA O (ASO)</b> <b>ANTY-STREPTOLYSIN O (ASO)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	138 IU/ml	110 – 166 IU/ml
<b>CZYNNIK REUMATOIDALNY (RF)</b> <b>RHEUMATOID FACTOR (RF)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	16.7 IU/ml	12.5 – 20.9 IU/ml

CORMAY IMMUNO-CONTROL I (L)

Nr kat. / Cat. No 4-288

Seria / Lot: 210-72

Data ważn. / Exp.: 2024.07

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
<b>BIALKO C-REAKTYWNE (CRP)</b> <b>C-REACTIVE PROTEIN (CRP)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	0.90 mg/dl	0.63 – 1.17 mg/dl
<b>ANTY-STREPTOLIZYNA O (ASO)</b> <b>ANTY-STREPTOLYSIN O (ASO)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	138 IU/ml	110 – 166 IU/ml
<b>CZYNNIK REUMATOIDALNY (RF)</b> <b>RHEUMATOID FACTOR (RF)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	16.7 IU/ml	12.5 – 20.9 IU/ml

CORMAY IMMUNO-CONTROL I (H)

Nr kat. / Cat. No 4-288

Seria / Lot: 210-70

Data ważn. / Exp.: 2024.07

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
<b>BIALKO C-REAKTYWNE (CRP)</b> <b>C-REACTIVE PROTEIN (CRP)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	2.50 mg/dl	2.13 – 2.88 mg/dl
<b>ANTY-STREPTOLIZYNA O (ASO)</b> <b>ANTY-STREPTOLYSIN O (ASO)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	301 IU/ml	256 – 346 IU/ml
<b>CZYNNIK REUMATOIDALNY (RF)</b> <b>RHEUMATOID FACTOR (RF)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	39.7 IU/ml	29.8 – 49.6 IU/ml

CORMAY IMMUNO-CONTROL I (H)

Nr kat. / Cat. No 4-288

Seria / Lot: 210-70

Data ważn. / Exp.: 2024.07

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
<b>BIALKO C-REAKTYWNE (CRP)</b> <b>C-REACTIVE PROTEIN (CRP)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	2.50 mg/dl	2.13 – 2.88 mg/dl
<b>ANTY-STREPTOLIZYNA O (ASO)</b> <b>ANTY-STREPTOLYSIN O (ASO)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	301 IU/ml	256 – 346 IU/ml
<b>CZYNNIK REUMATOIDALNY (RF)</b> <b>RHEUMATOID FACTOR (RF)</b>	TIA-wzmacniane lateksem Latex-enhanced TIA	39.7 IU/ml	29.8 – 49.6 IU/ml