



Myt-CAL

HEMATOLOGY CALIBRATOR

CALIBRATOR

LOT PLUS0424



2024-05-05

Instrument	Parameter	Calibration Values	Tolerance Limits
Orphée Mythic 18 II Generation	WBC x 10 ³ /μL	8.9	± 0.2
	RBC x 10 ⁶ /μL	4.52	± 0.10
	HGB g/dL	13.7	± 0.2
	HCT%	38.6	± 1.0
	PLT x 10 ³ /μL	233	± 12
Orphée Mythic 22	WBC x 10 ³ /μL	8.8	± 0.2
	RBC x 10 ⁶ /μL	4.44	± 0.10
	HGB g/dL	12.4	± 0.2
	HCT%	37.3	± 1.0
	PLT x 10 ³ /μL	245	± 12
Orphée Mythic 22 AL	WBC x 10 ³ /μL	8.8	± 0.2
	RBC x 10 ⁶ /μL	4.46	± 0.10
	HGB g/dL	12.5	± 0.2
	HCT%	37.2	± 1.0
	PLT x 10 ³ /μL	242	± 12
Orphée Mythic 60	WBC x 10 ³ /μL	8.9	± 0.2
	RBC x 10 ⁶ /μL	4.48	± 0.10
	HGB g/dL	12.8	± 0.2
	MCV fL	81.5	± 2.0
	PLT x 10 ³ /μL	250	± 12
	MPV fL	9.3	± 0.5
Orphée Mythic 70	WBC x 10 ³ /μL	9.1	± 0.2
	RBC x 10 ⁶ /μL	4.50	± 0.10
	HGB g/dL	12.7	± 0.2
	MCV fL	83.0	± 2.0
	PLT x 10 ³ /μL	258	± 12
	MPV fL	10.1	± 0.5

All brands and products are trademarks or registered trademarks of their respective companies.

